



# Glucometer Medisign MM1000

Code free, reliable and safe monitoring of blood glucose level in diabetes.

# **Glucometer Medisign MM1000**

## **Blood Glucose Monitoring System**

- · Continuously high measuring accuracy thanks to automatic double measurement
- Extra short measuring time (5 secs.)
- · Memory of 500 measuring values
- · Large illuminated display for easy reading of measurement values
- · Auto calibration, no coding required
- · With active sensor test strips
- CE marked
- Battery-operated\*

#### Accessories:

- 10 sensor test strips
- · Soft lancet device, soft lancets
- User manual
- Soft bag

## **Test strips**

- 50 sensor test strips
- Storage: 4°C 24°C
- Max. shelf-life 24 months
- CE marked

#### Available in two different units of measure:

- mml/L
- mg/dl



For more information, please contact us. We remain at your service!

action medeor e.V. German Medical Aid Organization

St. Töniser Str. 21 47918 Tönisvorst | Germany Phone: +49 2156 9788 100

info@medeor.de | www.medeor.de/en

